

Unique Glaucoma IOP Test Tonometer DIATON Presented at XXVII Congress of the ESCRS Barcelona, Spain 2009

Date: 09-19-2009 12:23 PM CET

Category: [Health & Medicine](#)

Press release from: [BiCOM Inc.](#)



Diaton tonometer is present at XXVII Congress of the European Society of Cataract and Refractive Surgeons ESCRS Meeting held in Barcelona, Spain. Tonometer presentation can be seen at booth E6 on the main floor.

Information on Diaton including procedure video can be found at TonometerDiaton.com

Unique Tonometer Diaton is presented at European Society of Cataract & Refractive Surgeons (ESCRS) meeting in the city of Barcelona, Spain. Transcleral and Trans-palpebral (through the Eyelid) Tonometry lets surgeons determine Intra-Ocular Pressure (IOP) shortly after surgery. The new Diaton Tonometer allows for quick non-contact screening of IOP through eyelid, independent of the cornea.

ESCRS is Europe's leading organization for cataract and refractive surgeons and is at the forefront of developments in anterior segment surgery. The society currently has 7,000 members from over 90 different countries worldwide. The event provided great exposure for Diaton Tonometer allowing for surgeons to fill in the gap.

Portability and simplicity make Tonometer Diaton ideal for a wide range of uses: mass screening of the population, at the patient's bedside, in geriatric homes, in children hospitals, for the military and for home use.

Diaton Tonometer is intended for use by Inpatient & Outpatient Clinics such as Hospitals, Emergency Rooms, Nursing & Elderly Homes, General & Specialty Practitioners as well as Ophthalmologists and Optometrists.

Tonometer Diaton is the perfect solution in the following cases when the use of other devices is problematic or impossible:

* mass prophylactic screening of patients;

- * IOP control during clinical observation of glaucoma patients;
 - * orthoclinostatical probe, as an additional test to diagnose glaucoma during selection of the adequate hypotensive therapy;
 - * ophthalmotone monitoring (even at night time);
 - * IOP measuring during contact correction (without need to remove contact lenses),
 - * IOP measuring in immobilized patients;
 - * IOP measuring in children.
- * on patients with the following conditions: chronic conjunctivitis, cornea pathology, including keratitis, keratone, cornea dimness, after penetrating keratoplastics, keratoprosthesis, laser refractive correction of the eyesight, high degree of ametropia, astigmatism;
- * on patients with medicinal allergies;
- * Lasik/ LASEK /PRK (recent clinical trials have proved that Diaton is the only device that can be used for IOP measurement right after these surgeries)

Besides determining the potential for glaucoma onset, Diaton tonometer vs other tonometers, Diaton can be used repeatedly during the day with no effect on the eye. The attendees of ESCRS are intrigued by the ease of use Diaton tonometer allows to measure IOP through the eyelid independent of the cornea.

BiCOM Inc. covers the full scope of customer service, which involves client training, follow ups, hotline support, sales support, etc. To get more information on the device or to make a purchase, visit the official website - www.TonometerDiaton.com, or call BiCOM at the toll free number 1-877-diatons (877.342.8667).

BiCOM Inc. , the global distributor of Tonometer Diaton is committed to preventing blindness worldwide through outreach, education and most importantly, mass screenings.

About BiCOM Inc.

BiCOM, with its 18 years of tradition and global clientele, has proven to be the enterprise of the highest level of professionalism, integrity and financial solvency. A unique team of engineers, medical, legal and business experts makes BiCOM Inc. the right place for global talent to find support and guidance. BiCOM Inc. sees its mission in bringing to market innovative products developed and manufactured worldwide.

more about Diaton tonometers www.TonometerDiaton.com

Press Contact –
Alina Lagoviyer / Media Relations
BiCOM Inc.
151 East Walnut Street
Long Beach, New York 11561
Phone: 1-877-342-8667
Fax: 1-888-260-0606
Contact@TonometerDiaton.com
www.TonometerDiaton.com

[You can find this press release here](#)